


¹³C-MQC Polarimeter



Specifically designed to interface with Oxford Instruments' Dynamic Nuclear Polariser (DNP) HyperSense[®], our new Polarimeter enables the user to routinely measure the level of polarisation in the liquid state sample after hyperpolarisation. This is particularly useful in the liquid state analysis of new compounds after hyperpolarisation, qualifying the polarisation performance of your HyperSense and in analyzing samples where a high-resolution spectrometer is not available.

A key feature of the Polarimeter is that it automatically calculates the level of sample hyperpolarisation, allowing rapid and consistent polarisation determination (calculation based on Pyruvic Acid).

The Polarimeter is a permanent magnet based, compact NMR spectrometer specially configured for calculating liquid state DNP.



The Polarimeter is a portable system with easy to use software, that can be located next to the DNP system or at a distance away.

After dissolution and transfer into the Polarimeter, the percentage liquid state polarisation is displayed. The Polarimeter has three calculation modes:

- Measurement of a reference sample for system calibration (normally ¹³C₂-DMSO)
- Measurement of polarisation enhancement of an aliquot (usually 0.5 to 4ml Pyruvic Acid solution) of a sample dissolved by HyperSense
- Measurement of all the sample dissolved by HyperSense to monitor the effective level of polarisation (¹³C₁ Pyruvic Acid)

Key enabling features of the Polarimeter are:

- Dedicated software which is easy to use by any operator
- Enhancements referenced to a standard so is independent of spectrometer performance
- Better SNR relative to the noisier thermal enhancement method at low NMR frequencies
- Larger 18mm probe so total sample volume is within the active coil area giving better accuracy

The Business of Science[®]

¹³C-MQC Polarimeter

Calculation overview

Monitoring the polarisation enhancement can be achieved by a number of methods, e.g. enhancement over thermal NMR spectrum, or polarisation spin population.

We have chosen to use spin polarisation as the benchmark because this can be calculated independently of the measuring spectrometer.

Example

A ¹³C₁ Pyruvic Acid (PA) sample is polarised whilst at the same time monitoring the solid state polarisation build-up in order to predict the maximum level of solid state polarisation.

After about 60 minutes the PA sample is dissolved into an 18mm NMR tube and the signal measured using the ¹³C-MQC Polarimeter. Then a sealed 4.0ml sample of ¹³C₂-DMSO is measured in the MQC and used as a comparison. From the integrals of these two samples the predicted liquid state polarisation at saturation is determined.




Figure 1: Hyperpolarised ¹³C NMR spectrum of PA (receiver gain 25) and ¹³C₂-DMSO (receiver gain 45) data obtained on the Polarimeter and processed.

¹³C-MQC Polarimeter Specifications

Use	Designed for use with the HyperSense DNP system
Electrical	Europe: 230V, 50Hz, USA: 110V, 60Hz, Japan: 100V+/- 10%, 50-60Hz
Magnet	0.55T permanent ¹³ C observe 18mm probe
Spectrometer	Single channel MQC spectrometer

RF Power amplifier

Single channel output	Nominal 300W at 6MHz 1% duty cycle pulsed at 100% power
Pulse length range	0.7µs- 1ms
Time resolution	0.1µs

Receiver

Gain	adjustable up to 72dB nominal in steps of 1dB from 12 to 72dB
Acquisition time	Maximum 80ms
ADC Resolution	16 bits

Dimensions

Magnet	W 362mm X D 344mm X H 342mm
Electronics unit	W 300mm X D 327mm X H 406mm
Weight	100Kg

Calibration sample

¹³C₂-DMSO 4.0ml

If you would like to know more about the Polarimeter or to discuss the application further please contact us.

This publication is the copyright of Oxford Instruments Molecular Biotools Ltd and provides outline information only, which (unless agreed by the company in writing) may not be used, applied or reproduced for any purpose or form part of any order or contract or regarded as the representation relating to the products or services concerned. Oxford Instruments' policy is one of continued improvement. The company reserves the right to alter, without notice the specification, design or conditions of supply of any product or service. HyperSense is a registered trademark of Oxford Instruments Molecular Biotools Ltd. and the Oxford Instruments logo is a trademark of Oxford Instruments plc or its subsidiaries.

HyperSense is not marketed for *in-vivo* use.

Polar02/08

Oxford Instruments Molecular Biotools Ltd.
Tubney Woods, Abingdon, Oxon OX13 5QX
Tel: +44 (0) 1865 393 200
Fax: +44 (0) 1865 393 333
E-mail: molecularbiotools@oxinst.co.uk
www.oxford-instruments.com

As part of
Oxford Instruments'
environmental policy
this brochure has been
printed on FSC paper

